KDB **GLIDE**®

Instructions For Use

CE
0459

Index

English (en)	1
Français (fr)	8
Deutsch (de)	16
Italiano (it)	24
Español (es)	32
Português (pt)	40
Čeština (cs)	48
Dansk (da)	56
Ελληνικά (el)	64
Suomi (fi)	72
Hrvatski (hr)	80
Latviski (lv)	88
Magyar (hu)	96
Nederlands (nl)	104
Polski (pl)	112
Română (ro)	120
Slovenščina (sl)	128
Svenska (sv)	136
Lietuvių k. (lt)	144
Norsk (no)	152
Íslenska (is)	160

A	Shaft
B	Indicator for Blade Orientation
C	Dual Blades
D	Pointed Tip
E	Heel
F	Ramp

 **CAUTION**

- A. Physician training is required before use of this instrument.
- B. Contents are guaranteed sterile provided that the package including the tray seal are not punctured or damaged. Do not use if the package is damaged or open, or if the seal appears damaged.  
- C. Read instructions before use.
- D. Do not use if the product is past its expiration date.
- E. When removing the instrument from its tray ensure that the tip does not touch any part of the tray. The blade tip may dull if it touches the tray surface. Do not touch the tip to any surfaces other than those of the eye as it may damage or dull the tip.
- F. Do not use instrument if it contains debris, is dull, bent, rusted or otherwise damaged.
- G. Precautions standard to ophthalmic surgery shall be observed when using the KDB GLIDE®, including customary aseptic technique.
- H. The instrument contains stainless steel. Some individuals may be sensitized or allergic to metals, avoid using in such cases.
- I. Used sharps are contaminated and can transmit diseases. Immediately after use, dispose of the instrument in an appropriately labeled puncture resistant biohazard sharps container.
- J. The KDB GLIDE® instrument is a single patient, single use disposable instrument. Do not reuse or re-sterilize. Re-sterilized or reused instruments pose risks of infection and dull blades that may result in damage to tissue.  

DESCRIPTION AND INTENDED USE

The KDB GLIDE® is a disposable ophthalmic knife, a manual surgical instrument used in ophthalmic surgical procedures to manually cut trabecular meshwork (TM) in pediatric and adult patients. The instrument is sterile and is intended for single patient, single use.

The KDB GLIDE® consists of a surgical grade stainless steel body and a plastic handle. The stainless steel body is comprised of a long shaft that allows for access to the TM, a pointed tip used to pierce TM, a ramp that lifts and stretches tissue as the device is advanced, and two blades that cut TM. An indicator on the plastic handle specifies the direction to advance the instrument.

INTENDED POPULATION

Adult and pediatric patients.

CONTRAINDICATIONS

Like most angle surgeries, proper visualization using a gonioscopy is key to manually cutting TM with the instrument. Do not attempt to use the instrument if corneal clarity is poor or if visualization of angle structures is not possible.

NOTICE: Any serious injury that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CAUTION: Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy or have poor visualization of angle structures. Improper visualization may result in misuse of the instrument and damage to eye structures adjacent to the TM.

UNDESIRABLE EFFECTS / SIDE EFFECTS

Adverse events that may be reasonably associated with the use of KDB GLIDE® in the eye include but are not limited to the following: elevation of intraocular pressure, hyphema, Descemet's membrane tear or detachment, iridodialysis, iris tear, iris injury, iritis, change to pupil shape, cyclodialysis cleft, corneal decompensation, corneal injury, corneal edema or opacification, perforation of sclera, as well as known

complications of general intraocular surgery including anterior chamber shallowing, hypotony, hypotony maculopathy, IOL dislocation, cataract formation, posterior capsule opacification, loss of vitreous, choroidal detachment, choroidal effusion, suprachoroidal hemorrhage, vitreous hemorrhage, posterior vitreous detachment, retinal detachment, retinal dialysis, retinal tear, endophthalmitis or other ocular infection.

The following Warnings, Precautions are measures to be taken in the event of malfunction of the device in its performance that may affect safety.

WARNINGS

- Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy.
- The KDB GLIDE® is intended to cut TM and cannot be used to make a clear corneal incision.

PRECAUTIONS

- Contents are provided sterile and intended for single patient, single use. Do not re-sterilize or re-use the device.
- Re-sterilized or reused instruments pose risks of infection and dull blades that may result in damage to tissue.
- Do not use if package, including the tray seal, are punctured or damaged. The KDB GLIDE® is intended to cut TM and cannot be used to make a clear corneal incision.
- Do not use past expiration date.
- When removing the instrument from its tray ensure that the tip does not touch any part of the tray. The blade tip may dull if it touches the tray surface. Do not touch the tip to any surfaces other than those of the eye as it may damage or dull the tip. Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy.
- Do not use instrument if it contains debris, is dull, bent, rusted, or otherwise damaged.
- Observe all usual precautions undertaken during intraocular surgery.
- Use aseptic technique and ensure KDB GLIDE® and field sterility as is customary during intraocular surgery.
- The instrument contains stainless steel. Some individuals may be sensitized or allergic to metals, avoid using in such cases.
- Used sharps are contaminated and can transmit diseases. Immediately after use, dispose of the instrument in an appropriately labeled puncture-resistant biohazard sharps container.

INSTRUCTIONS

1. Removing the Instrument from the Plastic Tray

- 1.1. **CAUTION:** See G.
- 1.2. Open the box and peel away the sealed Tyvek lid from the plastic container.
- 1.3. Remove the KDB GLIDE® from the plastic container taking care to lift the tip away from the tray first by gripping the KDB GLIDE® handle at the finger slots provided in the tray.
- 1.4. **CAUTION:** See E.
- 1.5. Inspect instrument. **CAUTION:** See F.

2. Considerations

- 2.1. Like other angle surgeries, use of a gonioscope over the patient's cornea under surgical microscope magnification is required for visualization of the TM and the instrument during use.
- 2.2. Positioning of the patient is essential to obtaining a proper view.
- 2.3. The KDB GLIDE® can be used stand-alone or combined with other procedures.
- 2.4. The KDB GLIDE® is inserted into the eye through a previously created clear corneal incision. Larger incisions allow for greater maneuverability.
CAUTION: The KDB GLIDE® is intended to cut TM and cannot be used to make a clear corneal incision.
- 2.5. The clear corneal incision should be approximately 180 degrees away from the TM intended to be cut.
- 2.6. The anterior chamber should be inflated with viscoelastic when using the KDB GLIDE®. **CAUTION:** Overinflating the chamber may make cutting of the TM difficult as the Canal of Schlemm may collapse.

4.

3. Using KDB GLIDE® under Gonioscopic View to cut TM.

- 3.1. Allowing some blood to reflux into the Canal of Schlemm, by lowering intraocular pressure, can be used to identify the TM. Insert the tip and shaft of KDB GLIDE® through the clear corneal incision and advance to the opposite side of the anterior chamber. **CAUTION:** Be careful to avoid touching the corneal endothelium and iris with the instrument as you are advancing or retracting the instrument in the anterior chamber. When inserting the tip, ensure that the long edge of the footplate is parallel to the incision (sharp tip is not perpendicular to incision slit). **(Figure a)**
- 3.2 Engage the TM across from the clear corneal incision with the pointed tip of the KDB GLIDE® and pierce through the TM and into the Canal of Schlemm while avoiding the anterior wall. **(Figure b)**
- 3.3 Advance KDB GLIDE® in the direction of the TM region that is intended to be cut by pivoting and advancing/retracting the shaft along the clear corneal incision while keeping the KDB GLIDE® footplate in the Canal of Schlemm. **(Figure c)** Avoid cutting surrounding tissue. Some blood reflux into the anterior chamber from the collector channels can be expected as the TM is removed. Viscoelastic or BSS can be used to improve visualization of the area being treated. Blood can be aspirated out of the anterior chamber. **CAUTION:** Ensure that the smooth heel of the tip slides forward with minimal pressure against the anterior wall of the Canal of Schlemm. If the heel is angled away from the back wall, the pointed tip may engage the anterior wall of the Canal of Schlemm and prevent smooth advancement of the KDB GLIDE® or result in sub-optimal engagement of the targeted TM tissue.
- 3.4 Detach the TM strip using KDB GLIDE® or microforceps. To detach the TM strip with the KDB GLIDE®, engage the TM across from the clear corneal incision with the pointed tip and pierce through the TM and into the Canal of Schlemm while avoiding the anterior wall of the Canal of Schlemm. Advance KDB GLIDE® in one direction until the desired endpoint of goniotomy is reached. Reverse the blade and advance towards the original area of treatment until a connection is made with the initial area of goniotomy.
- 3.5 Do not leave the detached TM strip in the anterior chamber. The detached TM tissue can be removed using KDB GLIDE® as the tissue may adhere to the blade, using ophthalmic forceps, or through aspiration from the clear corneal wound. Only aspirate fully detached TM strips.
- 3.6 **CAUTION:** If the angle of approach of the KDB GLIDE® instrument to the TM is too steep only one of the dual blades may incise the TM.
- 3.7 After cutting the TM, retract the KDB GLIDE® through the original clear corneal incision.

4. Disposal

- 4.1. **CAUTION:** See I.
- 4.2. **CAUTION:** See J.

HOW SUPPLIED

One KDB GLIDE® is supplied per box and comes packaged in a plastic tray with a sealed Tyvek lid. The tray secures and protects the delicate tip of the KDB GLIDE®. The KDB GLIDE® is a disposable, single patient, single use, sterile surgical instrument.

KDB GLIDE® is sterilized with gamma radiation.

STORAGE REQUIREMENTS

KDB GLIDE® should be stored at room temperature. **CAUTION:** Do not expose to water, sunlight, extreme temperatures, or high humidity.

EXPIRATION DATE

The KDB GLIDE® can be used up to 3 years after the date of manufacture as indicated by the expiration date on the label. Sterility is guaranteed until the expiration date if the packaging is not punctured or otherwise damaged. **CAUTION:** Do not use the KDB GLIDE® past the indicated expiration date.

Figure a




Figure b

6.

Figure c

USE OF SYMBOLS

The following symbols used in the KDB GLIDE® packaging are defined below.

	Manufacturer		Do not use if package is damaged
	Date of Manufacture YYYY-MM		Do not re-sterilize
	Distributor		Do not re-use
	"Use-by" or Expiration date		Consult Instructions for Use or Consult Electronic Instructions for Use https://www.newworldmedical.com/product-ifus/
QTY	Quantity of devices in package.		Caution
LOT	Lot Number		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
REF	Model/Catalogue Number	UDI	Unique Device Identifier
SN	Serial Number		Peel Direction
MD	Medical Device KDB GLIDE® Ophthalmic Instrument		Keep Dry
	Single Sterile Barrier System		Keep away from sunlight
	Single Sterile Barrier System with Protective Packaging Outside	EC REP	Authorized representative in the European Community
STERILE R	Sterilized using irradiation		

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Emergo Europe B.V.